

MEDICAL FITTINGS

A Division of Precision Medical, Inc.



Product Catalog
2011

TABLE OF CONTENTS

Flow Selector	2
Oxygen Timer	2
Ohmeda Style Quick Connects	3
Compact Ohmeda Style Check Units (Female Couplers)	4
Ohmeda Style Check Units (Female Couplers)	5
Chemetron Style Quick Connects	6
Chemetron Style Check Units (Female Couplers)	7
Puritan-Bennett Style Quick Connects	8
Puritan-Bennett Style Check Units (Female Couplers)	9
Schrader Style Quick Connects	10
Schrader Style Check Units (Female Couplers)	11
Oxequip Style Quick Connects	12
Medstar Quick Connects	13
"Y" Block Assemblies	14
DISS "T" and DISS "Y" Assemblies	14
T-Branch, Y-Block, Hex/Quad Block	15
Rail Blocks	15
Extension Manifolds	16
Pole Clamp	16
Ferrules	16
Medical Grade Hose	17
DISS Female	18
DISS Male	19
NPT Fittings	20-21
Miscellaneous	22-23
Condensation Traps	24
Neb Drives	25
International Fittings	26
Thread Comparsion Guide	28



FLOW SELECTOR and TIMOX

Max Pressure - 200PSI



Flow Selector

Flow Selector Features:

Give nebulizer treatments without disconnecting the humidifier.

Reduces changeover time to seconds.

No more loose tubing laying around.

Eliminates tubing nipple adaptors (it has 2 of its own).

Fits all flowmeters that have DISS male outlet fittings.

Lightweight.

Note: Flow is only delivered from one port at a time.



Typical Configuration

Flow Selector

PM1000



TIMOX
Oxygen Timer

TIMOX Features:

Automatic timing of oxygen flow as low as .25 lpm.

Unit is easily reset to zero with tip of pen or pencil.

Extremely accurate time measurement with time increments of .1 hours.

Large easy to read time display.

Up to 1999.9 recorded hours of oxygen flow.

Compact, lightweight and durable design.

Field-serviceable battery with an expected life of five years.

Two year warranty on unit excluding battery.

Timox Oxygen Timer 8050

Replacement Battery 8066

OHMEDA STYLE QUICK CONNECTS

(Available in ISO colors)

All Pipe Fittings are National Pipe Threads

Max Pressure - 200PSI

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD/Evac
 1/8 NPT Male	2101	2201	2301	2401	2601
 1/8 NPT Female	2102	2202	2302	2402	2602
 1/4 NPT Male	2103	2203	2303	2403	2603
 DISS Male	2105	2205	2305	2405	2605
 DISS Hex Nut	2106	2206	2306	2406	2606
 DISS Hand Tight	2108	2208	2308	2408	2608
 Hose Barb	1/4 Hose Barb 2111	1/4 Hose Barb 2211	1/4 Hose Barb 2311 5/16 Hose Barb 2312	1/4 Hose Barb 2411	1/4 Hose Barb 2611 5/16 Hose Barb 2615
 DISS Male with Check Valve	2114	2214	2314	2414	2614

COMPACT OHMEDA STYLE CHECK UNITS (Female Couplers)

(Available in ISO colors)

Max Pressure - 200PSI

All Pipe Fittings are National Pipe Threads

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide
 1/8 NPT Female	2160	2260	2360*	2460
 1/8 NPT Male	2161	2261	2361*	2461
 1/4 NPT Male	2163	2263	2363*	2463
 Hose Barb	1/4 Hose Barb 2164	1/4 Hose Barb 2264	1/4 Hose Barb 2364* 5/16 Hose Barb 2367*	1/4 Hose Barb 2464
 DISS Male	2165	2265	2365*	2465
 DISS Hex Nut	2166	2266	2366*	2466
 DISS Hand Tight	2168	2268	2368*	2468

*Vacuum Primary check only.

OHMEDA STYLE CHECK UNITS (Female Couplers)

(Available in ISO colors)

Max Pressure - 200PSI

All Pipe Fittings are National Pipe Threads

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD/Evac
 1/8 NPT Female	2150	2250	2350*	2450	2650
 1/8 NPT Male	2151	2251	2351*	2451	2651
 1/4 NPT Male	2153	2253	2353*	2453	2653
 Hose Barb	1/4 Hose Barb 2154	1/4 Hose Barb 2254	1/4 Hose Barb 2354* 5/16 Hose Barb 2357*	1/4 Hose Barb 2454	1/4 Hose Barb 2654 5/16 Hose Barb 8937
 DISS Male	2155	2255	2355*	2455	2655
 DISS Hex Nut	2156	2256	2356*	2456	2656
 DISS Hand Tight	2158	2258	2358*	2458	2658

*Vacuum Primary check only.

CHEMETRON STYLE QUICK CONNECTS

(Available in ISO colors)

All Pipe Fittings are National Pipe Threads

Max Pressure - 200PSI

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD/Evac
 1/8 NPT Male	3101	3201	3301	3401	3601
 1/8 NPT Female	3102	3202	3302	3402	3602
 1/4 NPT Male	3103	3203	3303	3403	3603
 DISS Male	3105	3205	3305	3405	3605
 DISS Hex Nut	3106	3206	3306	3406	3606
 DISS Hand Tight	3108	3208	3308	3408	3608
 Hose Barb	1/4 Hose Barb 3111	1/4 Hose Barb 3211	1/4 Hose Barb 3311 5/16 Hose Barb 3312	1/4 Hose Barb 3411	1/4 Hose Barb 3611 5/16 Hose Barb 3615
 DISS Male with Check Valve	3114	3214	3314	3414	3614

CHEMETRON STYLE CHECK UNITS (Female Couplers)

(ISO colors NOT available)

Max Pressure - 200PSI

All Pipe Fittings are National Pipe Threads

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD/Evac
 1/8 NPT Female	3150	3250	3350*	3450	3650
 1/8 NPT Male	3151	3251	3351*	3451	8940
 1/4 NPT Male	3153	3253	3353*	3453	3653
 Hose Barb	1/4 Hose Barb 3154	1/4 Hose Barb 3254	1/4 Hose Barb 3354* 5/16 Hose Barb 3357*	1/4 Hose Barb 3454	1/4 Hose Barb 3654 5/16 Hose Barb 8936
 DISS Male	3155	3255	3355*	3455	3655
 DISS Hex Nut	3156	3256	3356*	3456	3656
 DISS Hand Tight	3158	3258	3358*	3458	3658

*Vacuum Primary check only.

PURITAN-BENNETT STYLE QUICK CONNECTS

(Available in ISO colors)

Max Pressure - 200PSI

All Pipe Fittings are National Pipe Threads











	 Oxygen	 Air	 Vacuum	 Nitrous Oxide	 WAGD/Evac
 1/8 NPT Male	4101	4201	4301	4401	4601
 1/8 NPT Female	4102	4202	4302	4402	4602
 1/4 NPT Male	4103	4203	4303	4403	4603
 DISS Male with Knob	4105	4205	4305	4405	4605
 DISS Hex Nut	4106	4206	4306	4406	4606
 DISS Hand Tight	4108	4208	4308	4408	4608
 Hose Barb with Knob	1/4 Hose Barb 4111	1/4 Hose Barb 4211	1/4 Hose Barb 4311 5/16 Hose Barb 4312	1/4 Hose Barb 4411	1/4 Hose Barb 4611 5/16 Hose Barb 4615
 DISS Male with Knob and Check Valve	4114	4214	4314	4414	4614

PURITAN-BENNETT STYLE CHECK UNITS (Female Couplers)

(Available in ISO colors)

Max Pressure - 200PSI

All Pipe Fittings are National Pipe Threads


	 Oxygen	 Air	 Vacuum
 1/8 NPT Female	4150	4250	4350*
 1/8 NPT Male	4151	4251	4351*
 1/4 NPT Male	4153	4253	4353*
 Hose Barb	1/4 Hose Barb 4154	1/4 Hose Barb 4254	1/4 Hose Barb 4354* 5/16 Hose Barb 4357*
 DISS Male	4155	4255	4355*
 DISS Hex Nut	4156	4256	4356*
 DISS Hand Tight	4158	4258	4358*



*Vacuum Primary check only.

SCHRADER STYLE QUICK CONNECTS

All Pipe Fittings are National Pipe Threads

Max Pressure - 200PSI

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide
 1/8 NPT Male	6101	6201	6301	6401
 1/8 NPT Female	6102	6202	6302	6402
 1/4 NPT Male	6103	6203	6303	6403
 DISS Male with Knob	6105	6205	6305	6405
 DISS Hex Nut	6106	6206	6306	6406
 DISS Hand Tight	6108	6208	6308	6408
 Hose Barb	1/4 Hose Barb 6111	1/4 Hose Barb 6211	1/4 Hose Barb 6311	1/4 Hose Barb 6411
 DISS Male with Knob and Check Valve	6114	6214	6314	6414












Schrader Piston Identification	
	Hex Piston (standard)
	Swivel Piston (add "S" after part number)

SCHRADER STYLE CHECK UNITS (Female Couplers)

(Available in ISO colors)

Max Pressure - 200PSI

All Pipe Fittings are National Pipe Threads

	 Oxygen	 Air	 Vacuum	 Nitrous Oxide
 1/8 NPT Female	6150	6250	6350*	6450
 1/8 NPT Male	6151	6251	6351*	6451
 1/4 NPT Male	6153	6253	6353*	6453
 Hose Barb	1/4 Hose Barb 6154	1/4 Hose Barb 6254	1/4 Hose Barb 6354*	1/4 Hose Barb 6454
 DISS Male	6155	6255	6355*	6455
 DISS Hex Nut	6156	6256	6356*	6456
 DISS Hand Tight	6158	6258	6358*	6458












*Vacuum Primary check only.

OXEQUIP STYLE QUICK CONNECTS

(Available in ISO colors)

All Pipe Fittings are National Pipe Threads

Max Pressure - 200PSI












	 Oxygen	 Air	 Vacuum
 1/8 NPT Male	5101	5201	5301
 1/8 NPT Female	5102	5202	5302
 1/4 NPT Male	5103	5203	5303
 DISS Male with Knob	5105	5205	5305
 DISS Hex Nut	5106	5206	5306
 DISS Hand Tight	5108	5208	5308
 Hose Barb	1/4 Hose Barb 5111	1/4 Hose Barb 5211	1/4 Hose Barb 5311
 DISS Male with Knob and Check Valve	5114	5214	5314

MEDSTAR STYLE QUICK CONNECTS

(Available in ISO colors)

All Pipe Fittings are National Pipe Threads

Max Pressure - 200PSI

	 Oxygen	 Air	 Vacuum
 1/8 NPT Male	8101	8201	8301
 1/8 NPT Female	8102	8202	8302
 1/4 NPT Male	8103	8203	8303
 DISS Male	8105	8205	8305
 DISS Hex Nut	8106	8206	8306
 DISS Hand Tight	8108	8208	8308
 Hose Barb	1/4 Hose Barb 8111	1/4 Hose Barb 8211	1/4 Hose Barb 8311
 DISS Male with Check Valve	8114MS	8214	8314

"Y" BLOCK ASSEMBLIES, DISS "Y" and DISS "T" ASSEMBLIES

(Available in ISO colors)

Max Pressure - 200PSI




	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen
 Ohmeda Style	2152	2252	2352	2452	
 Compact Ohmeda Style	2162	2262	2362	2462	
 Chemetron Style	3152	3252	3352	3452	
 Puritan-Bennett Style	4152	4252	4352		
 Schrader Style	6152	6252	6352	6452	
 DISS "Y" Assemblies	7211	7212	7213	7214	7240
 DISS "T" Assemblies	7111	7112	7114		

T-BRANCH, Y-BLOCK, HEX/QUAD BLOCK and RAIL BLOCKS

Max Pressure - 200PSI

 <p>T-Branch</p>	<p>T-Branch, Standard 7100</p> <p>T-Branch, Offset 7100L/R</p>
 <p>Y-Block</p>	<p>1/8 NPT Female 7201</p> <p>1/4 NPT Female 7200</p>
 <p>Hex/Quad Block</p>	<p>1/4 NPT Female 9051</p>




Custom assemblies available upon request.

	Oxygen	Air	Vacuum																																				
 <p>Hex Block (1/8 NPT Female Inlet)</p>	7217	7218	7219																																				
 <p>Quad Block Assembly (1/8 NPT Female Inlet)</p>	7220	7221	7222																																				
 <p>Rail Block Assembly (DISS Male Inlet)</p>	<table border="1"> <thead> <tr> <th></th><th>Oxygen</th><th>Air</th><th>Vacuum</th></tr> </thead> <tbody> <tr> <td>Rail Block Assemblies*</td><td>7223</td><td>7224</td><td>7225</td></tr> <tr> <td colspan="4">*All Rail Block Assemblies includes Wall Bracket and Plate</td></tr> <tr> <td colspan="4">Replacement Parts</td></tr> <tr> <td>Mounting Bracket</td><td>7226</td><td></td><td></td></tr> <tr> <td>Mounting Plate</td><td>9056</td><td></td><td></td></tr> <tr> <td>Hex Block, 3 outlets</td><td>9054</td><td></td><td></td></tr> <tr> <td>Screws (2)</td><td>9058</td><td></td><td></td></tr> <tr> <td colspan="4">Custom assemblies available upon request.</td></tr> </tbody> </table>				Oxygen	Air	Vacuum	Rail Block Assemblies*	7223	7224	7225	*All Rail Block Assemblies includes Wall Bracket and Plate				Replacement Parts				Mounting Bracket	7226			Mounting Plate	9056			Hex Block, 3 outlets	9054			Screws (2)	9058			Custom assemblies available upon request.			
	Oxygen	Air	Vacuum																																				
Rail Block Assemblies*	7223	7224	7225																																				
*All Rail Block Assemblies includes Wall Bracket and Plate																																							
Replacement Parts																																							
Mounting Bracket	7226																																						
Mounting Plate	9056																																						
Hex Block, 3 outlets	9054																																						
Screws (2)	9058																																						
Custom assemblies available upon request.																																							

EXTENSION MANIFOLDS, POLE CLAMP and FERRULES


(Available in ISO colors)

Max Pressure - 200PSI

	Oxygen	Air
 Extension Manifold-2 Outlet	7230	7233
 Extension Manifold-3 Outlet	7231	7234
 Extension Manifold-4 Outlet	7232	7235
Extension Manifold Features: <ul style="list-style-type: none"> • Economical expansion of single wall outlet. • Permanently mounts to wall. • 2 foot conductive hose with 1/8 NPT Female inlet is standard. • Outlets are DISS Male with Check Valves. 		

 Pole Clamp	1" Pole 7216 1 1/2" Pole 7227
--	--

1/8 NPT Female holes (top and front) -
Comes with 1 plug on bottom

 Ferrule	1/4" Hose	
	1 Ferrule	9030
	25 Ferrules	9030-25
	50 Ferrules	9030-50
	5/16" Hose	
	1 Ferrule	503780
	25 Ferrules	503780-25
	50 Ferrules	503780-50

MEDICAL GRADE HOSE

Max Pressure - 200PSI



Kink resistant and light weight.

Reinforced with tough nylon braiding.

Color coded to conform to CGA-C9 standards.

Meets NFPA 99 requirements for conductivity
(*conductive hose only*).

FDA medical grade conductive hose.

Conveniently packaged in 250 foot rolls or as noted.

Maximum working pressure to 200 psi @ 70°F.

Non-Conductive Hose

Oxygen	Green	(1/4" ID-250' Roll)	9020-250
Air	Yellow	(1/4" ID-250' Roll)	9021-250
Vacuum	White	(1/4" ID-250' Roll)	9022-250
Nitrous Oxide	Blue	(1/4" ID-250' Roll)	9023-250
Nitrogen	Black	(1/4" ID-250' Roll)	9028-250
Carbon Dioxide	Gray	(1/4" ID-250' Roll)	9033-250
WAGD/Evac	Purple	(5/16" ID-250' Roll)	9038-250
	Clear	(1/4" ID-100' Roll)	504202-100

Conductive Hose

Oxygen	Green	(1/4" ID-250' Roll)	9120-250
Air	Yellow	(1/4" ID-250' Roll)	9121-250
Vacuum	White	(1/4" ID-250' Roll) (5/16" ID-250' Roll)	9122-250 9135-250
Nitrous Oxide	Blue	(1/4" ID-250' Roll)	9123-250
Nitrogen	Black	(1/4" ID-250' Roll)	9128-250
Carbon Dioxide	Gray	(1/4" ID-250' Roll)	9125-250
WAGD/Evac	Purple	(5/16" ID-250' Roll)	9124-250
	Black & White Striped	(1/4" ID-100' Roll)	9126-100








Custom Hose Assembly Ordering Information

Fitting	Length (Feet)	Type Hose	Fitting
—	—	—	—
<div style="border: 1px solid black; width: 100px; height: 30px; margin: 0 auto;"></div>	<div style="border: 1px solid black; width: 100px; height: 30px; margin: 0 auto;"></div>	<div style="border: 1px solid black; width: 100px; height: 30px; margin: 0 auto;"></div>	<div style="border: 1px solid black; width: 100px; height: 30px; margin: 0 auto;"></div>

DISS FEMALE

(Available in ISO colors)











Max Pressure - 200PSI

	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD/ Evac	Carbon Dioxide
 Hex Nut- 1/8 NPT with 1" Stem	0106						
 Hex Nut- 1/8 NPT with 1½" Stem	0107	0207	0307	0407	0507		504248
 Hex Nut- 1/4 NPT with 1½" Stem					0506	0606	
 Hex Nut- Hose Barb	1/4 Hose Barb 0115	1/4 Hose Barb 0216	1/4 Hose Barb 0316 5/16 Hose Barb 0317	1/4 Hose Barb 0416	1/4 Hose Barb 0516	1/4 Hose Barb 0616 5/16 Hose Barb 0617	1/4 Hose Barb 504154
 Hand Tight- 1/8 NPT with 1" Stem	0110						
 Hand Tight- 1/8 NPT with 1½" Stem	0111	0209	0309	0409	0509		
 Hand Tight- 1/4 NPT with 1½" Stem					0508	0608	
 Hand Tight- Hose Barb	1/4 Hose Barb 0118	1/4 Hose Barb 0218	1/4 Hose Barb 0318 5/16 Hose Barb 0319	1/4 Hose Barb 0418	1/4 Hose Barb 0518	1/4 Hose Barb 0618 5/16 Hose Barb 0619	1/4 Hose Barb 504287

DISS MALE

(Available in ISO colors)

Max Pressure - 200PSI

							
	Oxygen	Air	Vacuum	Nitrous Oxide	Nitrogen	WAGD/ Evac	Carbon Dioxide
DISS Male to 1/8 NPT Male	0701	0731	0751	0771	0791		
DISS Male with Check Valve to 1/8 NPT Male	0712	0742		0782	0797		
DISS Male to 1/4 NPT Male	0703	0733	0753	0773	0792	0794	8946
DISS Male with Check Valve to 1/4 NPT Male	0713	0743	0763	0783	0798	0799	8947
 1/4 Hose Barb	0709					0601 with Check 0602	
 1/4 Hose Barb with Knob	0710	0740	0760	0780	0795		
 1/4 Hose Barb with Knob & Check	0711	0741	0761	0781	0796		






NPT PIPE FITTINGS

Max Pressure - 200PSI

	1/8 NPT	1/4 NPT
 Pipe Tee	7030	7031
 Street Tee	7032	7033
 Female Tee	7034	7035
 Plug	7019	7017
 Drop Ear	Tee 7081 Elbow 7082	
 Cross	7036	7037
 45° Street Elbow	7056	7057
 90° Elbow	7058	7059
 90° Street Elbow	7063	7065
 Master Flow Valve	7064	7045

NPT PIPE FITTINGS

Max Pressure - 200PSI

 Hex Nipple	<table> <tr> <td>1/8 NPTM x 1/8 NPTM</td><td>7020</td></tr> <tr> <td>1/4 NPTM x 1/4 NPTM</td><td>7021</td></tr> <tr> <td>1/8 NPTM x 1/4 NPTM</td><td>7022</td></tr> </table>	1/8 NPTM x 1/8 NPTM	7020	1/4 NPTM x 1/4 NPTM	7021	1/8 NPTM x 1/4 NPTM	7022										
1/8 NPTM x 1/8 NPTM	7020																
1/4 NPTM x 1/4 NPTM	7021																
1/8 NPTM x 1/4 NPTM	7022																
 Connector	<table> <tr> <td>1/8 NPTF x 1/8 NPTM</td><td>1017</td></tr> <tr> <td>1/8 NPTF x 1/4 NPTM</td><td>7024</td></tr> <tr> <td>1/4 NPTF x 1/8 NPTM</td><td>7025</td></tr> <tr> <td>1/4 NPTF x 1/4 NPTM</td><td>7026</td></tr> <tr> <td>1/8 NPTF x 1/8 NPTM x 2.5" long</td><td>8955</td></tr> </table>	1/8 NPTF x 1/8 NPTM	1017	1/8 NPTF x 1/4 NPTM	7024	1/4 NPTF x 1/8 NPTM	7025	1/4 NPTF x 1/4 NPTM	7026	1/8 NPTF x 1/8 NPTM x 2.5" long	8955						
1/8 NPTF x 1/8 NPTM	1017																
1/8 NPTF x 1/4 NPTM	7024																
1/4 NPTF x 1/8 NPTM	7025																
1/4 NPTF x 1/4 NPTM	7026																
1/8 NPTF x 1/8 NPTM x 2.5" long	8955																
 Coupler	<table> <tr> <td>1/8 NPTF x 1/8 NPTF</td><td>7027</td></tr> <tr> <td>1/4 NPTF x 1/4 NPTF</td><td>7028</td></tr> <tr> <td>1/4 NPTF x 1/8 NPTF</td><td>7029</td></tr> </table>	1/8 NPTF x 1/8 NPTF	7027	1/4 NPTF x 1/4 NPTF	7028	1/4 NPTF x 1/8 NPTF	7029										
1/8 NPTF x 1/8 NPTF	7027																
1/4 NPTF x 1/4 NPTF	7028																
1/4 NPTF x 1/8 NPTF	7029																
 Hose Barb Connector	<table> <tr> <td>1/8 NPTM x 1/8" Dia.</td><td>7038</td></tr> <tr> <td>1/8 NPTM x 3/16" Dia.</td><td>7039</td></tr> <tr> <td>1/8 NPTM x 1/4" Dia.</td><td>7040</td></tr> <tr> <td>1/4 NPTM x 1/4" Dia.</td><td>7043</td></tr> <tr> <td>1/8 NPTM x 5/16" Dia.</td><td>8952</td></tr> <tr> <td>1/8 NPTF x 1/4" Dia.</td><td>7046</td></tr> <tr> <td>1/4 NPTF x 1/4" Dia.</td><td>7049</td></tr> <tr> <td>1/4 NPTF x 5/16" Dia.</td><td>8951</td></tr> </table>	1/8 NPTM x 1/8" Dia.	7038	1/8 NPTM x 3/16" Dia.	7039	1/8 NPTM x 1/4" Dia.	7040	1/4 NPTM x 1/4" Dia.	7043	1/8 NPTM x 5/16" Dia.	8952	1/8 NPTF x 1/4" Dia.	7046	1/4 NPTF x 1/4" Dia.	7049	1/4 NPTF x 5/16" Dia.	8951
1/8 NPTM x 1/8" Dia.	7038																
1/8 NPTM x 3/16" Dia.	7039																
1/8 NPTM x 1/4" Dia.	7040																
1/4 NPTM x 1/4" Dia.	7043																
1/8 NPTM x 5/16" Dia.	8952																
1/8 NPTF x 1/4" Dia.	7046																
1/4 NPTF x 1/4" Dia.	7049																
1/4 NPTF x 5/16" Dia.	8951																
 Pipe Nipple	<table> <tr> <td>1/8 NPTM x 1.5" Long</td><td>7050</td></tr> <tr> <td>1/8 NPTM x 2" Long</td><td>7051</td></tr> <tr> <td>1/8 NPTM x 2.5" Long</td><td>7052</td></tr> <tr> <td>1/4 NPTM x 1.5" Long</td><td>7053</td></tr> <tr> <td>1/4 NPTM x 2" Long</td><td>7054</td></tr> <tr> <td>1/4 NPTM x 2.5" Long</td><td>7055</td></tr> <tr> <td>1/8 NPTM x 3" Long</td><td>7062</td></tr> </table>	1/8 NPTM x 1.5" Long	7050	1/8 NPTM x 2" Long	7051	1/8 NPTM x 2.5" Long	7052	1/4 NPTM x 1.5" Long	7053	1/4 NPTM x 2" Long	7054	1/4 NPTM x 2.5" Long	7055	1/8 NPTM x 3" Long	7062		
1/8 NPTM x 1.5" Long	7050																
1/8 NPTM x 2" Long	7051																
1/8 NPTM x 2.5" Long	7052																
1/4 NPTM x 1.5" Long	7053																
1/4 NPTM x 2" Long	7054																
1/4 NPTM x 2.5" Long	7055																
1/8 NPTM x 3" Long	7062																
 Cap	<table> <tr> <td>1/8 NPTF</td><td>7018</td></tr> </table>	1/8 NPTF	7018														
1/8 NPTF	7018																
 Hose Splicer	<table> <tr> <td>1/4 Inside Diameter</td><td>7041</td></tr> </table>	1/4 Inside Diameter	7041														
1/4 Inside Diameter	7041																

MISCELLANEOUS







Max Pressure - 200PSI



 <p>0705 0706</p> <p>O₂ DISS 90° Elbow Male x Hex Nut</p>	<p>1" Stem 0705 1½" Stem 0706</p>
 <p>O₂ DISS Male x DISS Male</p>	0704
 <p>O₂ DISS Male x 1/8 NPT Female</p>	0702
 <p>O₂ DISS Female x O₂ DISS Female</p>	0104
 <p>O₂ Power Take-Off (1/8 NPT)</p>	0150
 <p>Air Power Take-Off (1/8 NPT)</p>	0151

MISCELLANEOUS

(Available in ISO colors)






Max Pressure - 200PSI

 <p>0154 0154C</p> <p>Hand Tight Tubing Nipple</p>	<p>Oxygen 0154</p> <p>Oxygen with Chain 0154C</p>
 <p>0158 0158C</p> <p>Wing Nut Tubing Nipple</p>	<p>Oxygen 0158</p> <p>Oxygen with Chain 0158C</p>
 <p>0159 0159C</p> <p>Hex Nut Tubing Nipple</p>	<p>Oxygen 0159</p> <p>Oxygen with chain 0159C</p>
 <p>Disposable Tubing Nipple</p>	<p>Clear 0119</p> <p>Yellow 0120</p> <p>Black 0121</p> <p>White 0122</p> <p>Green 0124</p>
 <p>9073 9090</p> <p>Tubing Nipple (for Vacuum Regulator)</p>	<p>Straight 9073</p> <p>90° 9090</p>
 <p>Chain Ring Assembly</p>	<p>6" Chain with Closed Ring 8004</p> <p>6" Chain with Open Ring 8005</p>

 <p>Crimping Tool</p>	9130
 <p>Liter Meter</p>	505652

CONDENSATION TRAPS

Max Pressure - 200PSI


	<p>Inlet DISS Female O₂ x Outlet DISS Male O₂ PM15-4</p>
	<p>Inlet DISS Female Air x Outlet DISS Male Air PM15-41</p>
	<p>Inlet DISS Female Air x Outlet DISS Male Air 90° Down PM15-43</p>
	<p>Inlet DISS Female O₂ 90° Up x Outlet DISS Male O₂ 90° Down PM15-44</p>
	<p>Inlet DISS Male Air 90° Down x Outlet DISS Female Air 90° Up PM15-45</p>

NEB DRIVE

(Available in ISO colors)











Max Pressure - 200PSI

		Oxygen	Air
 Ohmeda	6 to 7 lpm	8126	8226
	7 to 8 lpm	8127	8227
	8 to 9 lpm	8128	8228
	10 to 11 lpm	8129	8229
 Chemetron	6 to 7 lpm	8136	8236
	7 to 8 lpm	8137	8237
	8 to 9 lpm	8138	8238
	10 to 11 lpm	8139	8239
 DISS Hand Tight	6 to 7 lpm	8116	8216
	7 to 8 lpm	8117	8217
	8 to 9 lpm	8118	8218
	10 to 11 lpm	8119	8219
 Puritan-Bennett with Knob	6 to 7 lpm	8146	8246
	7 to 8 lpm	8147	8247
	8 to 9 lpm	8148	8248
	10 to 11 lpm	8149	8249
 Schrader	6 to 7 lpm	8166	8266
	7 to 8 lpm	8167	8267
	8 to 9 lpm	8168	8268
	10 to 11 lpm	8169	8269
 Oxequip with Knob	6 to 7 lpm	8156	8256
	7 to 8 lpm	8157	8257
	8 to 9 lpm	8158	8258
	10 to 11 lpm	8159	8259
 Medstar	6 to 7 lpm	8176	8276
	7 to 8 lpm	8177	8277
	8 to 9 lpm	8178	8278
	10 to 11 lpm	8179	8279

 Rescue Flow	Oxygen - 14 to 16 lpm		Oxygen - 24 to 26 lpm	
	DISS Hand Tight	8115	DISS Hand Tight	8114
	Ohmeda	8125	Ohmeda	8124
	Chemetron	8135	Chemetron	8134
	Puritan-Bennett	8145	Puritan-Bennett	8144
	Oxequip	8155	Oxequip	8154
	Schrader	8165	Schrader	8164
	Medstar	8175	Medstar	8174

INTERNATIONAL FITTINGS

Max Pressure - 200PSI

	Oxygen	Air	Vacuum	Nitrous Oxide
FRENCH (AFNOR)				
 Probe x 1/8 NPT Male	505979	505980	505981	
 Probe x 1/4 Hose Barb	505766	505767	505768	
EUROPEAN (NIST)				
 Nut/Nipple x 1/8 NPT Male	505752	505772	505773	
 Nut/Nipple x 1/4 Hose Barb	505774	505775	505776	
 Body Adaptor x 1/8 NPT Male	505777	505778	505779	
 Body Adaptor x 1/4 Hose Barb	505780	505781	505782	
SCANDINAVIAN (AGA)				
 Probe x 1/8 NPT Male	505763	505764	505765	
GERMAN (DIN)				
 Probe x 1/8 NPT Male	505753	505761	505762	
AUSTRALIAN (SIS)				
 Nut/Nipple x 1/8 NPT Male	505769	505770	505771	
BRITISH (BOC)				
 Probe x 1/8 NPT Male	6180	6184	6182	6183

SALES ORDER TERMS AND CONDITIONS

1. **TERMS OF PAYMENT:** Upon establishment of credit satisfactory to seller, terms of Net 30 Days are given unless otherwise stated. If credit is not established, payment must be received prior to shipment. Accounts not paid when due shall bear interest at the rate of one and one half percent (1-1/2%) per month until paid. Seller reserves the right at any time to suspend credit or to change credit terms provided herein when, in its sole opinion, the financial condition of Buyer so warrants. Failure to pay invoices at maturity date makes all subsequent invoices immediately due and payable, irrespective of terms, and Seller may withhold subsequent deliveries until the account has been settled in full. Applicant agrees to pay any collection costs incurred to collect any outstanding balance due, including reasonable attorney's fees. Buyer will be charged a fee for any returned checks.
2. **PRICES:** Prices quoted are firm for the duration of the original order. Prices are subject to review and adjustment in the event of Buyer initiated changes to the original order, whether written or oral, which materially affect part specifications, order quantity or schedule. Minimum order is \$50.00.
3. **MODIFICATIONS OR CANCELLATION:** Orders are not subject to modification or cancellation except with Seller's consent and subject to terms which indemnify Seller for all loss or damage.
4. **WARRANTY:** Precision Medical, Inc. warrants that its products will be free of defects in workmanship and/or material for the period specified in each product manual. Should any failure to conform to a product's warranty appear within the applicable period, Precision Medical, Inc. shall, upon written notification thereof and substantiation that the goods have been stored, installed, maintained and operated in accordance with Precision Medical, Inc.'s instructions and standard industry practice, and that no modifications, substitutions, or alterations have been made to the goods, correct such defect by suitable repair or replacement at its own expense. **ORAL STATEMENTS DO NOT CONSTITUTE WARRANTIES.** The representatives of Precision Medical, Inc. or any retailers are not authorized to make oral warranties about the merchandise described in this contract, and any such statements shall not be relied upon and are not part of the contract for sale. Thus, this writing is a final, complete and exclusive statement of the terms of that contract. **THIS WARRANTY IS EXCLUSIVE AND IS IN LIEU OF ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OTHER WARRANTY OF QUALITY, WHETHER EXPRESS OR IMPLIED.** Precision Medical, Inc. shall not under any circumstances be liable for special, incidental or consequential damages including but not limited to lost profits, lost sales, or injury to person or property. Correction of non-conformities as provided above shall constitute fulfillment of all liabilities of Precision Medical, Inc. whether based on contract, negligence, and strict tort or otherwise. Precision Medical, Inc. reserves the right to discontinue the manufacture of any product or change product materials, design, or specifications without notice.
5. **CLAIMS:** Any claims for erroneous changes, deficiencies or imperfections must be made within 30 days from the invoice date. Goods are not to be returned unless authorization has been secured from Seller.
6. **TITLE:** Unless otherwise agreed to in writing, all shipments are F.O.B. Seller's factory in Northampton, Pa and title to goods shall vest in Buyer upon delivery thereof to carrier.
7. **DELAYS:** Estimates of delivery are made by Seller in good faith but are not guaranteed. In the event of any delay in Seller's performance due to fires, flood, riot, strike or other difference with workman, accidents to machinery, acts of sabotage, failure of electrical power or other utilities, shortages of materials, labor or transportation, compliance with or other action taken to carry out the intent purpose of any law or regulation, or any cause beyond Seller's control, Seller shall have added time within which to perform this contract as may be reasonably necessary under the circumstance, and shall have the right to apportion its shipments among its customers in such a manner as it may consider equitable.
8. **TOOLING:** Specially designed tools and dies for which the Seller makes a separate charge remain the property of Seller and may not be removed by Buyer. Seller will maintain such tooling as long as needed to supply parts without further charge to Buyer and will not permit use by or for other than Buyer without Buyer's consent.
9. **PATENTS:** Buyer will indemnify, defend and save harmless Seller from liability for patent infringement on account of items manufactured to Buyer's specifications or sample.
10. **FAILURE OF SELLER TO ENFORCE RIGHTS:** Failure of Seller to enforce any of these terms and conditions or to exercise any right accruing from any default of Buyer, shall not affect or impair Seller's rights in case such default continues or in case of any subsequent default of Buyer.
11. **ERRORS:** Seller reserves the right to correct clerical or typographical errors without penalty.
12. **CONTROLLING PROVISIONS:** Seller shall not be bound by any terms on any Purchase Order or confirmation forms issued by Buyer with reference to the goods sold hereunder and which attempt to impose any conditions at variance with the terms stated herein. Seller's failure to object to provisions of such Buyer's form shall not be deemed a waiver or modification of the terms stated herein, which shall constitute the entire Agreement between parties. This Agreement shall be governed by and construed in accordance with the laws of the State of Pennsylvania.

Visit
www.precisionmedical.com


for
Suction Equipment
Compressors
Aspirators
Oxygen Regulators
Oxygen Conserving Regulators
Flowmeters

Thread Comparisons Guide:

(Actual Size shown)

Compare your old fittings to the threads to determine the NPT size you need.

(Measurements given are approximate.)





Medical Fittings

300 Held Drive

Northampton, PA 18067

Phone: 610-262-6020

Fax: 610-262-6080

Toll Free Phone: 800-331-2685

Toll Free Fax: 800-352-1240

Web Site: www.medicalfittings.com

ISO 9001 Certified