

CareFusion
Australia 316 Pty Limited
Unit 3/167
Prospect Highway
Seven Hills, NSW 2147
+61 2 9624 9000 tel
+61 2 9624 9034 fax

CareFusion
NZ 313 Ltd
14B George Bourke Drive
Mt. Wellington
Auckland
New Zealand
+64 9 270 2420 tel
0508 422 734 freephone
(within NZ)
+64 9 270 6285 fax

CareFusion
Unit 2 Oude Molen
Business Park
Oude Molen Road
Ndabeni
Cape Town 7405
South Africa
+27 21 510 7562 tel
+27 21 510 7567 fax

CareFusion
Hammarbacken 4B
SE-191 24 Sollentuna
Sweden
+46 8 544 43 200 tel
+46 8 544 43 225 fax

EC REP

CareFusion UK
The Crescent
Jays Close
Basingstoke
Hampshire RG22 4BS
United Kingdom
0800 917 8776 freephone
+44 (0)1256 330860 fax




CareFusion International Headquarters
A-One Business Centre
Zone d'activités Vers-la-Pièce no 10
1180 Rolle
Switzerland
+41 21 556 30 00 tel
+41 21 556 30 99 fax

CE
0086

© 2010 CareFusion Corporation or one of its subsidiaries. All rights reserved.
Alaris and Guardrails are registered trademarks of CareFusion Corporation
or one of its subsidiaries.
0000CF00304 Issue 1

Infusion Systems

Alaris® Products



carefusion.com



Infusion Systems

Advanced technology at the point of care

Take a close look at our product range and you can't fail to be impressed by the attention to detail. The Alaris® range of pumps has been specially developed to provide you with all the performance and functions your environment calls for.

Syringe Pumps

Alaris® CC Plus



Guardrails® Safety Software Ready	In-Line Pressure Monitoring	Drug Protocols	Alaris® Gateway Compatible	Theatre	Intensive Care	General Wards
-----------------------------------	-----------------------------	----------------	----------------------------	---------	----------------	---------------

✓	✓	✓	✓	✓	✓	✓
---	---	---	---	---	---	---

The top-of-the line Alaris® CC syringe pump with its variable pressure monitoring system is specially suited to help reduce some of the risks associated with undetected occlusions.

Alaris® GH Plus



✓		✓	✓		✓	✓
---	--	---	---	--	---	---

The Alaris® GH syringe pump with its advanced display provides comprehensive information at a glance, which is easy to access and change.

Alaris® PK



		✓	✓	✓	✓	
--	--	---	---	---	---	--

Specifically designed to meet the needs of the Anaesthetist, the Alaris® PK syringe pump contains software which includes Pharmacokinetic (PK) and Pharmacodynamic (PD) modelling for widely used anaesthesia drugs.

Alaris® TIVA



		✓	✓	✓	✓	
--	--	---	---	---	---	--

The Alaris® TIVA syringe pump combines ease-of-use with maximum control including the ability to pre-programme up to 100 protocols for Total Intravenous Anaesthesia (TIVA).

Alaris® Enteral



Enteral Milk Feeding	Neonatal
----------------------	----------

✓	✓
---	---

Designed to deliver Enteral therapy efficiently and safely.

Volumetric Pumps

Alaris® GP Plus



Guardrails® Safety Software Ready	In-Line Pressure Monitoring	Drug Protocols	Alaris® Gateway Compatible	Theatre	Intensive Care	General Wards
-----------------------------------	-----------------------------	----------------	----------------------------	---------	----------------	---------------

✓		✓	✓	✓	✓	✓
---	--	---	---	---	---	---

The Alaris® GP volumetric pump provides a flexible platform for IV therapy which is suitable across a wide range of applications.

Alaris® GW



Guardrails® Safety Software Ready	In-Line Pressure Monitoring	Drug Protocols	Alaris® Gateway Compatible	Theatre	Intensive Care	General Wards
-----------------------------------	-----------------------------	----------------	----------------------------	---------	----------------	---------------

			✓		✓	✓
--	--	--	---	--	---	---

Small, lightweight and easy to use, the Alaris® GW syringe pump is a versatile volumetric pump specially designed to fit within the Alaris® Patient Care System.

The Alaris® Editor

Improves consistency and standardisation.

Configuring pumps to specific care areas is easy and consistent with the Alaris® Editor. One PC tool configures the Alaris® CC and GH syringe pumps and the Alaris® GP volumetric pump.

- Centralised configuration control
- Up to 30 Care area profiles
- Up to 3000 drug protocols



‘Smart’ protection with Guardrails® safety software

Medication safety software is becoming a pre-requisite for infusion devices. However, many hospitals find the implementation a challenge. The **plus** software is designed to make the transition to full Guardrails® implementation easy and in line with your schedule.

You can use the same Alaris® Editor to create your drug library and care area profiles. When your Guardrails® data set is ready, upgrading the pump to Guardrails® is simple.



- Standardise**
Defining your infusion protocols
- Protect**
Your safety net
- Measure**
Powerful analysis of your infusion data
- Improve**
Drives Long term best practice of your evolving drug therapy