


# F&P Airvo 3


The Airvo 3™ is a respiratory support device that delivers warmed and humidified respiratory gases to spontaneously breathing patients through a variety of patient interfaces.

- A** Patient interface
- B** Pulse oximetry accessories
- C** High pressure oxygen (HPO) dual input manifold
- D** Oxygen-bottle holder
- E** Mobile pole stand
- F** Accessories basket
- G** Airvo 3 device
- H** Airspiral™ chamber
- I** Outlet elbow
- J** Battery module
- K** AirSpiral™ heated breathing tube
- L** Mobile pole stand handle




Airvo™ 3 system and compatible accessories

### Airvo 3 setup for NHF



**PT301xx\***  
Airvo 3



**900PT561**  
AirSpiral tube and chamber kit

Some products may not be available in your country. Please contact your local Fisher & Paykel Healthcare representative.

### Accessories



**900PT460D**  
**900PT460N**  
**900PT460S**  
High pressure oxygen (HPO) dual-input manifold (DISS, NIST, SIS)

**900PT462DN**  
HPO adapter (DISS to NIST)

**900PT474**  
Airvo 3 USB service cable

**900PT473**  
Airvo 3 data port adaptor



**900PT426**  
Storage basket



**900PT421**  
Mobile pole stand



**900PT428**  
C-clamp 25mm




**900PT427 (L)**  
Oxygen bottle holder




**900PT445**  
Mobile pole stand handle

### Disinfection




**900PT600**  
Disinfection kit



**900PT601**  
Disinfection filter



**900PT602**  
Cleaning sponge stick (20-pack)



**900PT603**  
Disinfection storage cover (20-pack)

### Spare Parts



**900PT448**  
Braked castor and fasteners



**900PT449**  
Steering castor and fasteners



**900PT450**  
Mobile pole stand hooks



**900PT930**  
Outlet elbow



**900PT957L**  
Battery module



**900PT933**  
Airvo 3 air filter

**900PT412xx\*** **900PT956**  
Power cord Power cord retainer

\* indicates a country code, if applicable.

Please note that the information on this material is summarized and provided for illustrative purposes only. For more information, please contact your local Fisher & Paykel Healthcare representative. The Airvo 3 and accessories are to be operated by, or under the supervision of, qualified personnel only. Refer to the manual and the instructions for use supplied with all accessories before using the device. F&P, Airvo, Optiflow and AirSpiral are trademarks of Fisher & Paykel Healthcare Limited. For patent information, see [www.fphcare.com/ip](http://www.fphcare.com/ip)